# Baumer Shoulder Prosthesis

Modular "MSP", "Neer" "TSS"

## DESCRIPTION:

A cemented prosthesis available in modular or non-modular conception, with anatomic design, which makes possible several types of "head x stem" assemblies, providing the surgeon with the best choice - according to the patient's characteristics and needs.

In the modular prosthesis, the humeral stem is available in four diameters (Ø 7 / Ø 9 / Ø 11 and Ø 13mm) with 125,0mm length. Humeral head in five different diameters (Ø 38/ Ø 41/ Ø 43/ Ø 45 and Ø 47mm). Concerned to non-modular prosthesis, the humeral stem is available in three diameters (Ø 8 / Ø 10 and Ø 12mm) in following different lengths: 120,0 / 130,0 and 140,0 mm, respectively. Humeral head in three diameters (Ø 34, Ø 40 and Ø 44mm). Glenoid components in two sizes (small and large). The stem and humeral heads are made of Titanium alloy Ti 6Al 4V EII which meets ASTM -136/ ISO 5832-3 specifications and Stainless Steel Cr Ni Mo ASTM F-138 / ISO 5832-1.

The glenoid component is fabricated of ultra-high -molecular-weight Polyethylene UHMW PE - ASTM F648 / ISO 5834-2.

Baumer shoulder prosthesis (when applicable) is coated with hydroxyapatite in the proximal region, by plasma spray method, with a thickness that provides a pore size of 50 to 150 microns.

## INDICATIONS:

- ..... Osteoarthritis
- ..... Posttraumatic osteoarthritis
- ..... Humeral head fracture sequel
- ..... Humeral head pseudoarthrosis
- ..... Revision of cemented arthroplasty

## HOW TO USE:

All the correct procedures and cares must be ensured when preparing the medullar canal as well as the glenoid region, to obtain perfect components positioning and fitting, avoiding radiolucency and the undesired micro-movements. Special instruments and tests to assist in this prosthesis surgical implantation are available. It is important that the instruments and tests used are those specifically designed for this system. Design variations, instrument sizes and similar tests may affect the critical measures required for an accurate implant.

Before starting the surgery, make sure that all the required instruments are available.

In order to protect the doctor and patient against possible future problems, it's recommended to write down on the patient's report the code and lot number of each implant used.

## CONTRAINDICATIONS:

- ..... All active infection and its sequels
- ..... All those known against the use of cemented prosthesis.

## PRECAUTIONS:

For safety and effectiveness in the use of this prosthesis, the surgeon must be acquainted with the procedures and recommendations for the implant (see information on how to use). To assure an adequate implant, only instruments and tests made by Baumer must be used. To maintain its sterility, the product should only be opened at the time of utilization and handled in sterile environments. In case of products supplied already sterilized, verify the sterilization validity. Do not use the implant if the package has been violated. In all the cases, traditional surgical practices must be followed in the postoperative period. The patient must be warned about the limitations involved in the total reconstruction of the joint and the recent history of this practice. Excessive physical activities and traumas, affecting the replaced joint, have implied the premature failure of the reconstruction by lack, fracture and/or wear of the implant. The patient must be warned about keeping the activities according to the situation and protecting the replaced joint against excessive stress. An implant should never be reused. Although it may be in good condition, the prior stress may have produced imperfections, that could reduce its useful life. The mixture of different materials may cause metallosis. In case of implant removal and forwarding for manufacturer analysis, it should be previously cleaned and disinfected by the user (surgeon/hospital/nurse). Therefore, the user should proceed with the methods mentioned below, or any other known and validated:

## Implantable Device/Disinfection:

Polymeric orthopedic components (UHMWPE): Cleaning method: Proteolytic enzyme solution, with ultrasonic treatment, or 3% stabilized hydrogen peroxide, or 1:100 solution sodium hypochlorite -Sterilizing or disinfection method: buffered, alkaline 2% solution of glutaraldehyde or ethylene oxide gas.

## Implantable Device/Disinfection:

Metallic orthopedic components:- Cleaning method: intense water rinse, 70% to 80% aqueous ethanol or isopropanol with subsequent ultrasonic treatment or proteolytic enzyme or 1:100 solution of sodium hypochlorite. Sterilizing or disinfection method: steam autoclave or ethylene oxide. To forward the removed implants to the manufacturer, packages must be used to retain the implant integrity. The package must contain information about the implant conditions, such as: used methods of cleaning and disinfection, as well as its identification.

## ADVERSE EFFECTS:

The implant of foreign materials in tissues may result in histological reactions involving several types of macrophages and fibroblasts. Cases of patients with reactions of sensitivity to metal after a total joint replacement have been rarely reported.

## STERILITY:

All the prosthesis components are supplied sterile. Either the metal components or those fabricated of UHMWPE are sterilized with ethylene oxide (EO). To resterilize unused components call our factory. All the instruments and tests must be sterilized before being used in a surgery. We recommend a steam autoclave. Never sterilize the UHMWPE components with vapor. Do not use the implant if the package has been violated. Check the sterilization validity date.

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